Appendix 2. ELISA working checklist

Dav	1
-----	---

and 630 nm.

1. Incubate the plate with 100 μL/well of **2 μg/mL EP1802Y Ab**. Seal and store at 4°C overnight. Day 2 1. Wash plate 3X with 300 µL/well of wash buffer and tap dry 2. Block by adding 250 µL/well of assay buffer to plate. Seal and incubate at RT for 1-3 hr on benchtop Start time: Sealed: → Stop time: 3. Prepare fresh standards from an aliquot of stock rPrP 4. Wash plate 3X with 300 μL/well of wash buffer and tap dry 5. Add 100 μL/well of rPrP standards, mouse QCs, and samples in duplicate. Seal and incubate at RT for 60-75 min. Start time: _____ → Stop time: Sealed: 6. Wash plate 3X with 300 µL/well of wash buffer and tap dry 7. Add 100 µL/well of **0.25 µg/mL biotin-8H4 Ab solution**. Seal and incubate at RT for 60-75 mins. Start time: → Stop time: Sealed: 8. Wash plate 3X with 300 µL/well of wash buffer and tap dry 9. Add 100 µL/well of **24.69ng/mL streptavidin-HRP solution**. Seal and incubate at RT for 30 mins. Start time: _____ → Stop time: 10. Wash plate 3X with 300 μL/well of wash buffer and tap dry. 11. Add 100 µL/well of RT **TMB**. Cover and incubate at RT on benchtop until Std. 1 (5ng/mL) reaches ~0.8 OD (pre-read at 605 nm) or 30 minutes max. Start time: _____ Covered: → Stop time: 12. Add 100 µL/well of RT Stop Solution. Mix well on plate reader briefly and read at 450 nm